
BIPLAST

**fr** Bande adhésive de contention  
**en** Adhesive compression bandage  
**de** Haftender Kompressionsverband  
**nl** Zelfklevend steunverband  
**it** Fascia adesiva di contenimento  
**es** Venda adhesiva de contención  
**pt** Banda adesiva de contenção  
**fi** Liimasiđe  
**sv** Självhäftande stödband  
**el** Αυτοκόλλητη ταινία σύγκράτησης  
**pl** Pas podtrzymujący samoprzylepny  
**lv** Liposa balsta saite  
**sk** Adhezívny fixačný obváz  
**hu** Öntapadós rugalmas fasíl  
**bg** Адезивен бинт за захранка  
**ro** Banda adeziva de compresie  
**ru** Бинт липкий фиксирующий  
**ko** 접착식 고정 봉대  
**ar** بند ملتصق لاصق

**fr** BANDE ADHÉSIVE DE CONTENTION**Description/Destination**

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Conserver cette notice.

**en****ADHESIVE COMPRESSION BANDAGE****Description/Destinazione**

Short-stretch adhesive compression bandage.

The device is available in different widths according to needs.

The device is intended to be used on the part of the body or the joint concerned by the patient's condition, and solely for the treatment of the listed indications.

**Composition**

Coton, polyester, adhésive · mass · (synthetic · elastomer · and · tackifying resin) · paper.

**Properties/Mode of action**

Adhesive compression bandage.

Aerated elastic fabric.

Multi-directional herringbone coating.

Short-stretch &lt; 70%.

Intermediate protective paper.

The device can be used for fixing dressings and splints.

**Indications**

Soft tissue or musculoskeletal system injuries.

Strapping: post-trauma support.

Fractures of the small bones (finger, foot, etc.).

Congenital or acquired bone deformity of the foot.

Vascular or lymphatic diseases requiring short-stretch or multicomponent bandage.

**Contraindications**

Do not use the product if the diagnosis has not been confirmed.

Do not use in the event of known allergy to any of the components.

Do not apply the product directly to an open wound without a dressing.

Do not use in case of oozing or major dermatological conditions of the affected limb.

Do not use in case of extra-anatomical bypass of the affected limb.

Do not use in case of severe peripheral neuropathy of the affected limb.

Do not use in case of septic thrombosis of the affected limb.

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**Precautions**

Verify the product's integrity before every use.

Do not use the product if it is damaged.

Strictly comply with your healthcare professional's prescription and recommendations for use.

If the skin is particularly sensitive, seek advice from a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

Remove the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Do not use after the expiry date indicated on the packaging.

**Undesirable side-effects**

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Any serious incidents related to the product should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

**Instructions for use/Application**

Retirez la pellicule de protection (Diagram ①).

Tear the paper to the desired length ①.

Pull each side to tear the entire length ②.

Peel off the paper on one side, apply to the limb, then peel off on the other side ③.

**Storage**

Store at room temperature, in the original packaging.

**Disposal**

Dispose of in accordance with local regulations.

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**THUASNE UK Ltd**
Thuasne SAS - SIREN/RCS Nanterre 542 091186