
THUASNE

Dynabelt

fr	Ceinture abdominale	7
en	Abdominal support	8
de	Bauchgurt	10
nl	Buikgordel	12
it	Cintura addominale	14
es	Faja abdominal	16
pt	Cinta abdominal	18
da	Mavebælte	20
fi	Tukivyö vatsalle	21
sv	Bukbälte	23
el	Ζώνη κοιλιακής χώρας	25
cs	Břišní pás	27
pl	Pas brzuszny	29
lv	Vēdera daļas josta	30
lt	Abdominalinis diržas	32
et	Kõhupiirkonna tugivöö	34
sl	Trebušni pas	35
sk	Brušný podporný pás	37
hu	Hasi támasztó öv	39
bg	Колан за коремната област	40
ro	Centură abdominală	42
ru	Бандаж для брюшного отдела	44
hr	Abdominalni pojas	46
zh	腹部支撐帶	48
ar	حزام للبطن	51

	1	66 - 80 cm
	2	81 - 90 cm
	3	91 - 100 cm
	4	101 - 120 cm*
	5	121 - 140 cm*
	6	141 - 160 cm*

fr	Tour de taille - Hauteurs : 16 cm/26 cm/33 cm *Pour la France uniquement, sur-mesure patient.
en	Waist measurement - Heights: 16 cm/26 cm/33 cm *For France only, made to measure patient.
de	Taillenumfang - Höhen: 16 cm/26 cm/33 cm *Nur für Frankreich, maßgeschneidert.
nl	Tailleomvang - Hoogtes: 16 cm/26 cm/33 cm *Enkel voor Frankrijk, op maat van de patiënt.
it	Circonferenza vita - Altezze: 16 cm/26 cm/33 cm *Solo per la Francia, su misura per il paziente.
es	Contorno de cintura - Alturas: 16 cm/26 cm/33 cm *Solo para Francia, a medida del paciente.
pt	Perímetro de cintura - Alturas: 16 cm/26 cm/33 cm *Apenas para a França, à medida do paciente.
da	Taljemål - Højde: 16 cm/26 cm/33 cm *Kun for Frankrig, efter patientens mål.
fi	Vyötärönympäry - Korkeudet: 16 cm / 26 cm / 33 cm *Ainoastaan Ranskassa, potilaan mittojen mukaan.
sv	Midjemått - Höjd: 16 cm/26 cm/33 cm *Endast för Frankrike, patientspecifik.
el	Περίμετρος μέσης - Ύψη: 16 εκ./26 εκ./33 εκ. * Μόνο για Γαλλία, κατά παραγγελία ασθενούς.
cs	Obvod pasu - Výšky: 16 cm / 26 cm / 33 cm *Výhradně pro Francii, na míru pacienta.
pl	Obwód w talii - Wysokości: 16 cm/26 cm/33 cm *Tylko we Francji produkt jest wykonywany na miarę dla pacjenta.
lv	Vidukļa apkārtmērs - augstumi: 16 cm/26 cm/33 cm *Tikai attiecībā uz Franciju - pacientam pielāgots izstrādājums.
lt	Juosmens apimtis - Aukštis: 16 cm / 26 cm / 33 cm *Tik Prancūzijai, pagal paciento užsakymą.
et	Vööümberrõõt - Kõrgused: 16 cm/26 cm/33 cm * Ainult Prantsusmaal, patsiendile kohandatud.
sl	Obseg pasu - Višina: 16 cm/26 cm/33 cm *Samo za Francijo, po meri pacienta.
sk	Obvod pása - Výšky: 16 cm/26 cm/33 cm *Len pre Francúzsko, vyrobený na mieru pacienta.
hu	Derékbőség - Magasság: 16 cm/26 cm/33 cm *Csak Franciaországban, a páciens egyéni méretére szabva.
bg	Обиколка на талията - Височини: 16 cm/26 cm/33 cm *Само за Франция, изработен по мярка за пациента.
ro	Circumferința taliei - Înălțimi: 16 cm/26 cm/33 cm *Doar pentru Franța, făcut pe măsura pacientului.
ru	Обхват талии - Высота: 16 см/26 см/33 см * Только для Франции, изготавливается по меркам пациента.
hr	Opseg struka - Visine: 16 cm/26 cm/33 cm *Samo u Francuskoj, po mjeri bolesnika.
zh	腰围-高度 : 16 cm/26 cm/33 cm *仅限于法国，为患者定制。
ar	محيط الخصر - الارتفاعات: 16 سم / 26 سم / 33 سم *موجه لفرنسا فقط، حسب مقياس المريض.

Compression

Soutien

Compression

Support

Kompression

Stütze

Compressie

Ondersteuning

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Sostegno

Compresión

Sujeción

Compressão

Suporte

Kompression

Støtte

Kompressio

Tuki

Kompression

Stöd

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ABDOMINAL SUPPORT

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table*.

Composition

Aerated and breathable Combitex fabric. **(1)**

Closure with self-fastening surface and thumb-loops. **(2)**

Flexible and removable dorsal stays. **(3)**

Textile components: coated elastane - polyamide - polyester.
Non-textile components: brushed tempered steel - polyamide.

Properties/Mode of action

The device ensures strong support and adaptable compression to the abdomen and the lumbar area thanks to the elastic fabric et the flexible stays.

Indications

Abdominal muscle weakness.

Postoperative conditions: incisional hernia, umbilical hernia or epigastric hernia, laparotomy, laparoscopy, Laparo-endoscopic single site surgery, cesarean section.

Postpartum conditions.

Bruising to the rib cage.

Rib fracture.

Minor upper and lower back pain.

Contraindications

Do not apply the product in direct contact with broken skin.

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of a hiatus hernia.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

Do not use in the event of Pelvic floor deficiency.

Surgeons usually do not recommend any kind of compression for a period of time after cardiac or thoracic surgery for people at risks.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

It is recommended that a healthcare professional adjust the stays to the shape of the patient's back before the belt is worn for the first time.

Do not wear the product in a medical imaging machine.

Do not wear the product while sleeping.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional. For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support without excessive compression.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

The use of a medical device by a child, should be done under the supervision of an adult or a healthcare professional.

Undesirable side-effects

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

In postoperative conditions, the restrictive effect of the device could alter respiratory function and decrease lung capacity and tidal volume.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Position the belt around the waist, using the thumb-loop. ④

The belt must be centred on the spine: the dorsal stays should be positioned on either side of the spine.

The bottom of the belt must stop at the top of the crease between the buttocks.

The inside label must be positioned upwards.

Hold the left end flat using your left hand.

Pull the thumb-loop simultaneously and evenly. ⑤

Fasten the belt, right side over left side, adjusting to the desired compression. ⑥

If the self-fastening strap is properly positioned it does not wrinkle.

The level of compression can be adjusted as required throughout the day.

Care

Close the self-fastening tabs before washing. Machine washable at 60 °C. Do not machine wash the device more than 10 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat. Can be sterilised in an autoclave at 120°C.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

*Not applicable to made-to-measure medical devices, where measurements will be taken individually for each patient.

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Zakopčajte pojas na desnoj strani preko lijeve strane prilagođavajući željenu razinu zatezanja. ©
Thuasne SAS - SIREN/RCS Nanterre 542 091 186 -