THUASNE®

Genuimmo 0°

fr	Attelle d'immobilisation du genou à 0°	4
en	Knee immobilisation splint at 0°	4
de	Schiene zur Ruhigstellung des Knies in einem 0°-Winkel	4
nl	Immobilisatiebrace voor de knie op 0°	5
it	Tutore di immobilizzazione per il ginocchio a 0°	5
es	Férula de inmovilización de la rodilla a 0°	6
pt	Tala para imobilização do joelho em 0°	6
da	Immobiliserende knæskinne ved 0°	6
fi	Immobilisoiva polvituki, taivutus 0°	7
sv	Immobiliseringssskena för knät 0°	7
el	Νάρθηκας ακινητοποίησης γονάτου σε κάμψη 0°	8
cs	Imobilizační ortéza kolene 0°	8
pl	Orteza do unieruchomienia kolana pod kątem 0°	8
lv	Celgala imobilizācijas saite 0°	9
lt	Kelio imobilizavimo 0° padėtyje įtvaras	9
et	Põlve fikseerimise lahas 0° juures	10
sl	Opornica za imobilizacijo kolena s kotom 0°	10
sk	Imobilizačná kolenná ortéza pri 0° uhle	10
hu	Oldalsínes térdrögzítő 0°-os szög	11
bg	Шина за имобилизация на коляно на 0°	11
ro	Atelă pentru imobilizarea genunchiului la 0°	12
ru	Шина для имобилизации коленного сустава под углом 0°	12
hr	Udlaga za imobilizaciju koljena, 0°	13
zh	0°膝关节固定夹板	13
ar	جيرة لتثبيت الركبة على زاوية صفر درجة	13

		< 175 cm	> 175 cm
	S	35 - 45 cm	x
	M	45 - 52 cm	
	L	53 - 62 cm	
	XL	63 - 75 cm	

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - 使用指南 - إرشادات الاستعمال


© Thuasne - 6173501 (2023-07)
© Studio Caterin


THUASNE
120, rue Marius Aujan
92300 Levallois-Perret
France

www.thuasne.com








			< 175 cm	> 175 cm
	S	35 - 45 cm	 50 cm	x
	M	45 - 52 cm		 60 cm
	L	53 - 62 cm		
	XL	63 - 75 cm		

fr	Taille patient	Circonférence cuisse	Hauteur
en	Patient height	Thigh circumference	Height
de	Größe des Patienten	Oberschenkelumfang	Höhe
nl	Maat van de patiënt	Dijomtrek	Hoogte
it	Altezza del paziente	Circonferenza coscia	Altezza
es	Altura del paciente	Circunferencia del muslo	Altura
pt	Tamanho do paciente	Circunferência da coxa	Altura
da	Patientens højde	Lårets omkreds	Højde
fi	Potilaan koko	Reiden ympärysmitta	Korkeus
sv	Patientens storlek	Omkrets lår	Höjd
el	Μέγεθος ασθενούς	Περιφέρεια μηρού	Ύψος
cs	Výška pacienta	Obvod stehna	Výška
pl	Rozmiar pacjenta	Obwód uda	Wysokość
lv	Pacienta augums	Augšstilba apkārtmērs	Augstums
lt	Paciento dydis	Šlaunies apimtis	Aukštis
et	Patsiendi suurus	Reie ümbermõõt	Kõrgus
sl	Velikost pacienta	Obseg stegna	Višina
sk	Veľkosť pacienta	Obvod stehna	Výška
hu	A beteg mérete	Combkörméret	Magasság
bg	Размер на пациента	Бедрена обиколка	Височина
ro	Talia pacientului	Circumferință coapsă	Înălțime
ru	Размер пациента	Окружность бедра	Высота
hr	Veličina pacijenta	Opseg bedra	Visina
zh	患者身高	大腿圍	高度
ar	مقاس المريض	محيط الفخذ	الارتفاع

Immobilisation
Immobilisation
Ruhigstellung
Immobilisatie
Immobilizzazione
Inmovilización
Imobilizaçào
Støtte
Immobilisoi
Immobilisering
Ακινητοποίηση
Znehybnění
Unieruchomienie
Imobilizācija
Imobilizavimas
Liikumatuks muutmine
Imobilizacija
Imobilizácia
Rögzítés
Обездвижване
Imobilizare
Имобилизация
Imobilizacija
固定
التثبيت

fr

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en
KNEE IMMOBILISATION SPLINTS AT 0° AND 20°

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

The device is composed of:

- 3 removable, radiolucent, rigid stays,
- 4 self-fastening closure and support straps.

Option of opening the side of the splint while still keeping the knee immobilised (see application).

Knee splint with fine fabric and well-ventilated thanks to large openings.

No seams against the skin.

Bilateral product.

Available in 4 sizes for the short version and 3 sizes for the long version.

Composition

Textile components: PVC - polyamide - cotton - polyurethane - polyester - natural latex - elastane.

Rigid components: aluminium - silicone - acetal.

Properties/Mode of action

Knee joint stabilisation via the rigid stays.

Total immobilisation at 0° or 20°.

Adjustable closure and support via self-fastening straps.

Indications

Post-operative immobilisation.

Post-traumatic immobilisation.

Knee pain.

Contraindications

Do not use the product if the diagnosis has not been confirmed. Do not use in the event of known allergy to any of the components. Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Do not use in the event of unstable fractures, fractures of the lower part of the tibia or fibula, or fractures of the upper part of the femur.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

It is recommended that the venous thromboembolic risk be assessed before any immobilisation is indicated. Follow the advice of a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

In case of history of venous or lymphatic disorders, consult a healthcare professional.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Prolonged immobilisation can sometimes cause muscle weakness. Do not wear the device when driving a vehicle.

Contains natural latex.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

Choose the splint height, then the size.

Measurements should be taken at the widest part of the thigh. It is recommended that a healthcare professional supervises the first application.

It is recommended to wear the device directly on the skin, unless contraindicated.

- Open the splint and lay it flat open.
- Place the opened splint under the leg, the wider edge towards the top of the leg.

N.B.: The top of the splint must not impede flexion of the leg at the groin.

• Shaping the stays: if needed, it is possible to shape the stays in order to avoid potential pressure points or the splint slipping.

Apply a slight pressure in the desired direction.

Repeat the process on the other stays if needed.

If subsequent shaping of the stays is needed, consult a healthcare professional.

• Tighten the straps in the order described.

Tightening of ① and ② by passing them under the lateral uprights of the splint.

Fastening of the panels ③ and ④.

Tightening of ⑤ and ⑥ by passing through the buckles.

Once all the straps are fastened, readjust the tightening.

• Splint adaptability:

The straps located just above and below the knee are height-adjustable ⑦, meaning that they can be adapted to fit the patient's shape and condition (presence of scar following surgery, wound, etc.).

The knee splint can be opened on the outside ⑧, particularly in the event of surgery.

To do this, simply unfasten the hook located in the middle of the splint and pass it under the upright before fastening it at the back of the knee splint.

If needed, reduce the strap length by cutting the edges with scissors.

Care/Maintenance

Product can be washed in accordance with the instructions shown on this leaflet and on the label. Machine washable at 30°C (delicate programme). Close the self-fastening tabs before washing. Remove stays before washing. Stays can be washed with a moist cloth.

Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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Thuasne UK Ltd