
**THUASNE®**

# HipLoc Evo

de Gebrauchsanweisung - Hüftgelenkorthese  
en Instructions for use - Hip joint orthosis  
fr Consignes d'utilisation - Orthèse de hanche  
it Istruzioni per l'uso - Ortesi per l'articolazione dell'anca  
es Instrucciones de utilización - Órtesis de cadera  
pt Instruções de uso - Ortótese para a articulação da anca  
nl Gebruiksaanwijzing - Heuporthese  
da Brugsanvisning - Hofledsortose  
sv Bruksanvisning - Höftledsortos  
fi Käyttöohje - Lonkkanivelen tukiside  
el Οδηγίες χρήσης - Νάρθηκας άρθρωσης ισχίου  
cs Návod k použití - Ortéza kyčelního kloubu  
pl Sposób użytkowania - Orteza stawu biodrowego  
et Kasutusjuhend - Puusaliigese ortoos  
lv Lietošanas instrukcija - Gūžas locītavas ortoze  
lt Naudojimo instrukcija - Klubo sąnario įtvaras  
sk Návod na použitie - Ortéza bedrového kĺbu  
sl Navodila za uporabo - Ortoza za kolčni sklep  
hu Használati utasítás - Csípőízület ortézis  
ro Instrucțiuni de utilizare - Orteză articulație șold  
bg Инструкции за употреба - Ортеза за тазобедрена става  
hr Upute za uporabu - Ortoza za zglob kuka

<sup>1</sup>Hip bandage · Ceinture à la taille · Tutore dell'anca · Cincha de la cadera · Bandagem da anca · Heupbandage · Hoftebandage · Höftbandage · Lantiovyö · Κηδεμόνας ισχίου · Kyčelní bandáž · Bandaż biodrowy · Puusabandaaz · Gūžas bandāža · Klubo tvarsčiai · Bedrová bandáž · Bandaža kolka · Csíróbandázs · Bandaj pentru șold · Бандаж за бедро · Bandaža kuka

<sup>2</sup>Pocket for pelvic shell · Poche pour la coque pelvienne · Tasca per l'elemento di sostegno bacino · Bolsillo para la tablilla pélvica · Bolsa para soporte da bacia · Zak voor de bekkenschaal · Lomme til bækkenskål · Påse för bäckenskål · Lantiokuoren tasku · Θήκη για το πυελικό κέλυφος · Obal na kyčelní podpěru · Kieszeń dla łuski biodrowej · Seljatoe tasku · Kabata iegurņa ieliktnim · Dubens apvalkalo jđeklas · Obal bedrovej podpery · Žep za medenično čašico · Táška medencehéj szárnára · Buzunar pentru cochilie șold · Джеб за тазовия компонент · Džep za bedrenu ljusku

<sup>3</sup>Rear section · Partie dos · Parte posteriore · Parte trasera · Parte das costas · Ruggpand · Rygdæl · Ryggdæl · Selkäosa · Μέρος της πλάτης · Zádová část · Element tylny · Seljaosa · Muguras daļa · Nugarinė dalis · Chrbtová časť · Hrbtnei del · Hátrész · Partea posterioră · Задна част · Leđni dio

<sup>4</sup>Thigh pad · Coussinet de la cuisse · Pelotta per la coscia · Almohadilla de la muslera · Almofada da coxa · Dijbeenpelotte · Lårpelotte · Lårpelott · Reisipehmuste · Μαξιλαράκι μηρού · Stehenní pelota · Pelota udowa · Reiepelott · Augšstilba spilventiņš · Paminkštinta šlaunų pagalvėlė · Stehenná pelota · Stegenska blazinica · Combárnázat · Pelotă pentru coarsă · Подложка за бедрото · Bedrena pelota

<sup>5</sup>Splint system · Atelle · Sistema di stecche · Sistema de guías · Sistema de talas · Spalk · Skinnesystem · Skensystem · Lasta · Σύστημα με μπανέλες · Systém dlah · System szyn · Lahasesüsteem · Šinu sistēma · Skydelių sistema · Systém dláh · Sistem tirnic · Sínrendszer · Sistem cu șine · Шинна система · Sustav udлга

<sup>6</sup>Cover for splint system · Housse pour l'attelle · Copertura per il sistema di stecche · Forro para el sistema de guías · Revestimiento para sistema de talas · Hoes voor de spalk · Betræk til skinnesystem · Överdrag för skensystem · Lastan suojus · Κάλυμμα για το σύστημα με μπανέλες · Poťah na systém dlah · Pokrowiec dla systemu szyn · Lahasesüsteemi kate · Šinu sistēmas pārvalks · Skydelių sistemės įvalkalas · Poťah na dlahy · Preveleka za sistem tirnic · Huzat sínrendszer szárnára · Îneltoare sistem cu șine · Кальф за шинната система · Presvlaka za sustav udлга

<sup>7</sup>Thigh shell · Coque de la cuisse · Guscio tutore coscia · Muslera · Suporte da coxa · Dijbeenschaal · Lårskål · Lårskål · Reisiikuori · Κέλυφος μηρού · Stehenní podpěra · Łuska udowa · Reietugi · Augšstilba ieliktnis · Šlaunies apvalkalas · Stehenná podpera · Stegenska čašica · Combhéj · Cochilie pentru coarsă · Компонент за бедрото · Bedrena ljuska

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Oberschenkelpelotte: ABS - Acryl-Butadien-Styrol

Bandagensystem: PA - Polyamid

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## **Mobilising hip joint orthosis with adjustable extension and flexion**

### **Purpose of use**

HipLoc Evo is an orthosis for mobilisation of the hip joint with adjustable extension and flexion. The orthosis may only be used for hip treatment.

### **Mode of action**

The HipLoc Evo orthosis is used as a luxation protection after endoprothetic hip joint replacement. By means of a controlled approval of movements, the patient is induced to avoid any movements tending to cause luxation.

### **Indications**

- Instability of the hip joint
- Stabilisation of the hip joint after total endoprothetic replacement with loose muscle control and a tendency to luxation
- Stabilisation of the hip joint after repositioning a luxated hip joint following endoprothetic replacement

### **Contraindications**

In case of the following medical conditions, consult a doctor before using the orthosis:

- Injuries and skin disorders, serious scars with swelling, reddening skin and increased temperature in the orthosis contact area
- Severe varicose illnesses (varicose veins) with venous insufficiencies
- Severe varicose illnesses (varicose veins) with venous insufficiencies
- Sensitivity disorders in the legs (e.g. in case of diabetes mellitus)
- Circulation disorders

Never allow the product to come into contact with injured skin.

### **Side effects**

No side effects are known if properly used.

### **Description of the orthosis**

The orthosis features a splint and bandage system in which the splint system of the orthosis comprises a pelvic shell on one side, the splint system, a thigh pad and a thigh shell.

The splint system here connects the one-sided pelvic shell with the thigh shell; and in the hip joint area of which a joint is integrated with an adjustable extension and flexion unit. The splint system is adjustable and may be individually adjusted to the patient's specific leg length. If necessary, the splint system can be provided with a cover which may be variably cut to size.

The pelvic shell (thermoplastically remouldable) is located on the side and dorsally encloses the pelvis up to sacrum height. The joint should ideally be positioned to match the position of the anatomical transverse hip joint axis. Flexion can be limited to 90° and 60°, and extension to 0°.

With its movable suspension system, the thigh shell adjusts to the leg movements, thus ensuring that the shell always sits flush on the leg.

A cushioned thigh pad is fixed to the thigh splint on the body side; the pad is intended to exert lateral pressure on the thigh muscles and improves guidance of the femur. The pressure exerted should, however, not be uncomfortable for the patient. The pad's position may be adjusted in increments of 2 cm.

The bandage system consists of non-elastic and elastic materials and may be individually adjusted for waists of between 80 – 140 cm. A wide hook and loop fastener in the front area additionally allows the individual closure of the fastener to any anatomical characteristic.

### **Fitting the orthosis**

#### **Hip bandage:**

1. Shorten the bandage to the measured waist circumference. The extension integrated in the back section can be variably shortened with a pair of scissors. After shortening, simply fasten the end of the bandage again to the rear side element – by the hook and loop system. **Tip:** In order to improve the fit and provide a waist cut of the belt, the extension can be fastened at an angle to the rear section using the hook and loop system.

2. Open the hook and loop fastener of the pocket provided on the outside of the bandage and insert into it the pelvic brace which is already connected with the splint system. Thereafter, fixate the hook and loop fastener again. (Should the splint system not yet be connected with the pelvic brace, please connect both parts with each other prior to insertion into the pocket.)

#### **Splint system:**

3. In order to adjust the flexion and extension limit, simply loosen the pins in the pelvic shell and reposition them accordingly.

4. Adjust the splint system by pushing relative to one another the two splints sitting on top of each other and, at the right length, fix them into place with two screws. The thigh shell should here be positioned above the knee. The closure of the thigh shell can be variably shortened and adjusted to any leg circumference.

**Caution:** Should the thigh shell not yet be connected with the splint system, both elements must first be screwed together. To this end, open the pocket with its hook and loop fastener on the

outer side of the thigh shell and connect it with the lower end of the splint. Close the pocket subsequently again with the hook and loop fastener. If necessary, the aluminium splints can be put crosswise and adjusted to match the patient's body contours.

5. The thigh pad should now be repositioned in its height. First remove the pad cover and then loosen the two screws and washers (between pad and splint); then position and fasten the pad at the desired height. After that, pull the pad cover over the pad again.

6. If necessary, a cover can be pulled over the adjusted splint system. In order to pull over the cover, simply remove the thigh shell and pad from the splint system and push the cover from below with its perforated side towards the leg via the splints. The overhang on the lower part of the cover can be fixed to the fastener point provided at the lower end of the splint. If necessary, the upper part of the cover may be shortened in its length. After that, fasten again the thigh shell and pad to the splint system.

### Use

The period of use of the orthosis thus prepared is to be set after consultation with the treating physician.

### Special note

All screws must be locked in place after fitting (e.g. using Loctite or similar).

### Precautions

Follow the advice of the professional who prescribed or supplied the product. Consult your professional if discomfort occurs. Store at room temperature, preferably in its original packaging. For hygiene and performance reasons, do not reuse the product for another patient. Any serious incident that occurs in connection with the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient are/is domiciled.

### Care

- Always wash the orthosis separately and do not soak it
- Do not use any cleaning agents, conditioners or aggressive products (containing chlorine)
- Remove the splint system and the thigh shell prior to washing and close the hook-and-loop fastener systems
- The splint system and the thigh shell may be cleaned with a damp cloth
- Squeeze out the water thoroughly
- Dry away from any heat sources

### Material composition

Pelvic and thigh shell: PE - polyethylene

Guidance system: aluminium

Thigh pad: ABS - acrylonitrile butadiene styrene

Bandage system: PA – polyamide

Please keep this instruction leaflet in a safe place.

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Thuasne Scandinavia AB · Box 95 · 131 07 Nacka
**Thuasne UK Ltd**
**Thuasne Deutschland GmbH**