THUASNE®

Ligaflex® Immo

fr	Attelle d'immobilisation du genou - 3 volets.....	3
en	Knee immobilisation splint - 3 panels.....	3
de	Schiene zur Ruhigstellung des Knies - 3 Teile.....	3
nl	Immobilisatiebrace voor de knie - 3 flappen.....	4
it	Tutore di immobilizzazione per il ginocchio - 3 pannelli.....	4
es	Férula de inmovilización de la rodilla - 3 volets.....	4
pt	Tala para imobilização do joelho - 3 painéis.....	5
da	Immobiliserende knæskinne - 3 fløje.....	5
fi	Immobilisoiva polvituki - 3 läppää.....	5
sv	Immobiliseringsskena för knät - 3 delar.....	6
el	Νάρθηκας ακινητοποίησης γονάτου - 3 πτερύγες.....	6
cs	Imobilizační ortéza kolene - 3dílné provedení.....	7
pl	Orteza do unieruchomienia kolana - 3 składowe elementy.....	7
lv	Celgala imobilizācijas saite - 3 daļas.....	7
lt	Kelio imobilizavimo įtvartis - 3 dalių.....	8
et	Põlve fikseerimise lahas - 3-osaline.....	8
sl	Opornica za imobilizacijo kolena - 3 deli.....	8
sk	Imobilizačná kolenná ortéza - 3-dielna.....	9
hu	Oldalsínes térdrögzítő - 3 részes.....	9
bg	Шина за имобилизация на коляно - 3 части.....	9
ro	Atelă pentru imobilizarea genunchiului - 3 componente.....	10
ru	Шина для иммобилизации коленного сустава - 3 вставки.....	10
hr	Udlaga za imobilizaciju koljena - 3 krilca.....	10
zh	膝关节固定夹板 - 3片式.....	11
ar	جبيرة لتثبيت الركبة - 3 مفصلات.....	11

	1	50 - 60 cm
	2	60 - 70 cm
	3	70 - 80 cm
	4	> 80 cm

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung -
Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso -
Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet -
Råd vid användning - Οδηγίες χρήσης - Pokyny k používání - Szczegóły
użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos -
Kasutusjuhised - Navodila za uporabo - Pokyny na používanie -
Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare -
Инструкция по эксплуатации - Upute za uporabu - 使用指南 -
ارشادات الاستعمال

© Thuasne - 6172901 (2023-07)
© Studio Caterin

THUASNE
120, rue Marius AUFAN
92300 Levallois-Perret
France



www.thuasne.com

	1	50 - 60 cm
	2	60 - 70 cm
	3	70 - 80 cm
	4	> 80 cm

fr	Hauteur entrejambe - sol
en	Inside leg measurement
de	Innenahthöhe - Boden
nl	Hoogte kruis - grond
it	Altezza cavallo
es	Altura de entrepierna - suelo
pt	Altura entre as pernas-solo
da	Skridthøjde - gulv
fi	Korkeus haara-alue - lattia
sv	Höjd innerben - golv
el	Ύψος από τη βουβωνική χώρα έως το έδαφος
cs	Vzdálenost mezi rozkrokem a zemi
pl	Wysokość od kroku do podłoża
lv	Kājas iekšējais augstums - zeme
lt	Vidinės siūlės ilgis nuo grindų
et	Kõrgus hargmikust tallani
sl	Notranja dolžina noge - tla
sk	Dĺžka od podložky po rozkrok
hu	Föld -lábközép magasság
bg	Височина от горната част на бедрото до пода
ro	Înălțimea coapsă - sol
ru	Длина ноги по внутренней стороне
hr	Visina u preponama - od tla
zh	裤裆高度-地面
ar	ارتفاع الساق من الداخل - الأرض

Immobilisation
Immobilisation
Ruhigstellung
Immobilisatie
Immobilizzazione
Inmovilización
Imobilização
Støtte
Immobilisoi
Immobilisering
Активтоποίηση
Znehybnění
Unieruchomienie
Imobilizācija
Imobilizavimas
Liikumatuks muutmine
Imobilizacija
Imobilizácia
Rögzítés
Обездвижване
Imobilizare
Имобилизация
Imobilizacija
固定
التثبيت

fr

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KNEE IMMOBILISATION SPLINT - 3 PANELS

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

The device is composed of :

- 4 removable, radiolucent, rigid stays,
- 3 panels, adjustable to various leg circumferences,
- 7 self-fastening closure and support straps,
- foam padding for patient comfort.

Bilateral product.

Available in 4 sizes.

Composition

Textile components : polyester - polyurethane - cotton - polypropylene - polyamide.

Rigid components : aluminium - polyamide.

Properties/Mode of action

Knee joint stabilisation via the rigid stays.

Total immobilisation at 0°.

Adjustable closure and support via self-fastening straps.

Indications

Post-operative immobilisation.

Post-traumatic immobilisation.

Knee pain.

Contraindications

Do not use the product if the diagnosis has not been confirmed. Do not use in the event of known allergy to any of the components.

Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Do not use in the event of unstable fractures, fractures of the lower part of the tibia or fibula, or fractures of the upper part of the femur.

Precautions

Verify the product's integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

It is recommended that the venous thromboembolic risk be assessed before any immobilisation is indicated. Follow the advice of a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

In case of history of venous or lymphatic disorders, consult a healthcare professional.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product's performance, remove it and consult a healthcare professional.

Prolonged immobilisation can sometimes cause muscle weakness. Do not wear the device when driving a vehicle.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters, etc.) or wounds of various degrees of severity.

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended that a healthcare professional supervises the first application.

It is recommended to wear the device directly on the skin, unless contraindicated.

- Open the splint and lay it flat open **1**.
- Place the opened splint under the leg **2**, the wider edge towards the top of the leg. The patella must be centred in the cut-away area. N.B. : The top of the splint must not impede flexion of he leg at the groin.
- Shaping the stays: if needed, it is possible to shape the stays in order to avoid potential pressure points or the splint slipping. Apply a slight pressure in the desired direction. Repeat the process on the other stays if needed. If subsequent shaping of the stays is needed, consult a healthcare professional.
- Position the side panels in order to fully cover the leg (panels are adjustable to the thigh and calf circumference) **3**.
- Engage the self-fastening straps in the buckles. Tighten the straps in the order described (A, B, C, D, E, F, G) **4**. Once all the straps are fastened, readjust the tightening. If needed, reduce the strap length by cutting the edges with scissors.

Care/Maintenance

Product can be washed in accordance with the instructions shown on this leaflet and on the label. Hand wash. Close the self-fastening tabs before washing. Remove stays before washing.

Stays can be washed with a moist cloth. Do not use detergents, fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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Thuasne UK Ltd