



RETAILER IDENTIFICATION

Patient's surname:

Case No. for renewal

Patient's first name:

Gender : M F Child

1st treatment

Patient's height:

Date : Quantity:

Customer code:

I authorize my health care professional to collect my data and to communicate them to Thuasne company as part of the processing of my made to measure medical device in accordance with Law No 78-17 of 6 January 1978, and European Regulation No 2016/679/ EU of 27 April 2016, I have rights including in particular the rights of access, rectification, portability and deletion of my data. I can exercise these rights by contacting the health care professional to who I ordered my medical device.

Patient signature

If possible, please enclose photos of the limb to be fitted.

Please draw in the contours of the garment on the diagram and cross unnecessary measures.

RIGHT ARM LEFT ARM

Fill in one form for each side

Indications

- Major burn (by default)
- EDS (Ehlers-Danlos Syndrome)

Fabrics

- Cicatrex Nature
- Cicatrex Filifine
- Cicatrex AirSkin

	WHITE	BEIGE	BLACK
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-	<input type="checkbox"/>	<input type="checkbox"/>

Options

Sleeve fastening

- Without anti-slip
- Elastic plain braid 3 cm
- Anti-slip with silicone dots 3 cm
- Silicone men anti-slip 5 cm
- Silicone women anti-slip 5 cm
- Up to the armpit
- Shoulder cover
- Shoulder attachment (bra loop)

Opening

- Type:
- Zip
 - Velcro

Dimensions: cm

Position:

- Join side
- Back of the hand side
- Palm of the hand side

Other

- Compression pad

Dimensions (l x h): cm

Position: to be drawn in

Comments

.....

.....

Heights	Measures in cm
V Length from the wrist to the middle of the widest part of the forearm	
W Length from the wrist to the elbow	
X Length from the wrist to the middle of the widest part of the upper arm	
Y Length from the wrist to the armpit	
H Length from the armpit (g) to the acromion process (tip of the shoulder)	
Circumferences	
C Wrist circumference	
D Circumference at the middle of the widest part of the forearm	
E Elbow circumference	
F Circumference at the middle of the widest part of the upper arm	
G Arm circumference at the armpit	
I Chest circumference from the acromion process and going under the opposite armpit	

Circumferences in cm Length in cm

