



- ORDER (by default)
- QUOTATION
- RENEWAL



RETAILER IDENTIFICATION

Customer code:

Patient's surname: .....

Patient's first name: .....

Gender :  M  F  Child

Patient's height: .....

Case No. for renewal


1st treatment

Date : ..... Quantity: .....

I authorize my health care professional to collect my data and to communicate them to Thuasne company as part of the processing of my made to measure medical device in accordance with Law No 78-17 of 6 January 1978, and European Regulation No 2016/679/ EU of 27 April 2016, I have rights including in particular the rights of access, rectification, portability and deletion of my data. I can exercise these rights by contacting the health care professional to who I ordered my medical device.

Patient signature

**If possible, please enclose photos of the limb to be fitted.**  
**Please draw in the contours of the garment on the diagram and cross unnecessary measures.**  
**Please also fill in torso form No. 1b.**



○ circumference in cm    □ length in cm    ◡ half-circumference in cm

35 and 36 : For the body only  
 Measurement taken at an angle for a brief-cut body  
 Measurement taken straight for a short-cut body

**Please contact your regular Thuasne distributor**



RETAILER IDENTIFICATION

**Patient's surname:** .....  
 Patient's first name: .....  
 Gender :  M  F  Child  
 Patient's height: .....

Case No. for renewal

1st treatment

Date : ..... Quantity: .....

**Customer code:**

I authorize my health care professional to collect my data and to communicate with me to Thuasne company as part of the processing of my made to measure medical device in accordance with Law No 78-17 of 6 January 1978, and European Regulation No 2016/679/ EU of 27 April 2016, I have rights including in particular the rights of access, rectification, portability and deletion of my data. I can exercise these rights by contacting the health care professional to who I ordered my medical device.

**Patient signature**

**If possible, please enclose photos of the limb to be fitted.  
 Please draw in the contours of the garment on the diagram and cross unnecessary measures.  
 Please also fill in torso form No. 1a.**

### Indications

- Major burn (by default)  
 EDS (Ehlers-Danlos Syndrome)

### Models

- Shoulder pad  
 Body  
 Vest

### Fabrics

- Cicatrex Nature  
 Cicatrex Filifine  
 Cicatrex AirSkin

	WHITE	BEIGE	BLACK
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	-	<input type="checkbox"/>	<input type="checkbox"/>

### Options

#### Center opening

- Zip  
 Front position  
 Back position  
 Hook to help fasten the vest (by default)

#### Armpit compression

- Right  Yes  No\*  
 Left  Yes  No\*

\* Soft fabric by default

### Collar

- Round neck  
 "V" neck  
 Front collar cut height (from the sternum):  cm  
 Back collar cut height (from the 7th cervical vertebra):  cm  
 Stand-up collar  
 Circumference at top:  cm  
 Circumference at bottom:  cm  
 Front height:  cm  
 Back height:  cm  
 Top of sternum/neck height:  cm

### Cups

- With compression  
 Without compression  
 Depth:  A  B  C  D  E  F

### Velcro fastening (for panty)

- Vest over panty  
 Panty over vest

### Additional opening\*\*

- Position: (to be drawn in)  
 Zip  Velcro Length:  cm

### Soft fabric\*\*

- Elbow  
 Other .....

### Other\*\*

- Compression pad  
 Position: (to be drawn in)  
 Dimensions:  cm  
 Belt

\*\* Extras, not included in the cost

**RIGHT SIDE**

13 Shoulder length from the acromion process to neck circumference No. 1

14 Armhole circumference

16 Top of sleeve

15 Arm circumference at the armpit

18 Length to the middle of the upper arm

17 Upper arm circumference (widest part)

20 Length to the elbow

19 Elbow circumference

22 Length to the middle of the forearm

21 Forearm circumference (widest part)

24 Length to the wrist

23 Wrist circumference

36 Thigh circumference\*

1 Neck circumference from the 7th cervical vertebra to the top of the sternum

2 Chest width, breathing in (armpit to armpit)

4 Top of sternum - armpit line length

3 Chest perimeter, breathing in

6 Top of sternum - under the bust length

5 Circumference under the bust, breathing in

8 Top of sternum - waist length

7 Waist circumference

10 Top of sternum - bottom of vest length (for a body, length to the crotch)

9 Hip circumference

37 Back length (7th cervical vertebra - bottom of vest). For a body, length to the crotch, measure with the tape measure flat against the skin

39 Circumference tip of the breast breathing in

38 Length top of sternum - tip of the breast

a\*\* Length neck base - tip of the breast

b\*\* Length neck base - under the bust

**LEFT SIDE**

12 Armhole circumference

11 Shoulder length from the acromion process to neck circumference No. 1

25 Arm circumference at the armpit

26 Top of sleeve

27 Upper arm circumference (widest part)

28 Length to the middle of the upper arm

29 Elbow circumference

30 Length to the elbow

31 Forearm circumference (widest part)

32 Length to the middle of the forearm

33 Wrist circumference

34 Length to the wrist

35 Thigh circumference\*

\* For the body only  
 \*\* To fill out only from D cup

**Please contact your regular Thuasne distributor**