



THUASNE LYMPHOLOGY STOCKINGS

POSITION OF THE ANATOMICAL MARKERS

- G - is located at the gluteal fold (for the tights)
- F - is equidistant between G and E
- E - corresponds to the middle of the patella
- D - is located on the head of the fibula (3 - 4 cm below the patella)
- C - corresponds to the widest point of the calf
- B1 - corresponds to the start of the calf muscles
- B - to be measured at the narrowest point of the ankle
- Y - goes through the instep and the tuberosity of the calcaneus
- A - goes through the head of the metatarsal bones (at the widest point)



All markers are positioned on the external face of the leg



Metatarsal head circumference


Instep circumference

Circumference of the ankle (narrowest point)

Circumference of the thigh (at gluteal fold)

Measuring heights (position of measuring tape on the external face of the leg)

Buttocks circumference



Crotch : Distance at the front (K1T) / Distance at the back (K2T)

Hip circumference

Waist circumference (or indicate the circumference at the desired height)



RETAILER IDENTIFICATION

Patient's surname:

Case No. for renewal

Patient's first name:

Gender : M F Child

1st treatment

Patient's height:

Date : Quantity:

Customer code:

I authorize my health care professional to collect my data and to communicate them to Thuasne company as part of the processing of my made to measure medical device in accordance with Law No 78-17 of 6 January 1978, and European Regulation No 2016/679/ EU of 27 April 2016, I have rights including in particular the rights of access, rectification, portability and deletion of my data. I can exercise these rights by contacting the health care professional to who I ordered my medical device.

Patient signature

**If possible, please enclose photos of the limb to be fitted.
Please draw in the contours of the garment on the diagram and cross unnecessary measures.**

Models

- Below-knee
- Thigh-high
- Half-tight
- Tights

LEFT RIGHT

Compression

- Class 2 (15 - 20 mmHg)
- Class 3 (20 - 36 mmHg)
- Class 4 (> 36 mmHg)

LEFT RIGHT BRIEF*

* The compression level of the briefs must be identical to one of the two legs.

Below-knee options

- Anti-slip with silicone dots 3 cm
- Anti-slip with silicone dots 5 cm
- Without anti-slip (ribbed edge)

LEFT RIGHT

Thigh-high options

- Without anti-slip
- Silicone anti-slip for women 5 cm
- Silicone anti-slip for men 5 cm
- Anti-slip with silicone dots 5 cm

LEFT RIGHT

Half-tight options

- Short part side

Short part measure (ℓGF):

LEFT RIGHT

 cm

General options

Toe

- Open toe
- Closed toe

LEFT RIGHT

Compression pad

- Instep compression pad
- Medial malleolus compression pad
- Lateral malleolus compression pad

Dimensions (ℓ x h):


 cm

Foam type

- Soft 7 mm
- Hard 7 mm

Comments

.....



Please contact your regular Thuasne distributor