
MOBIDERM SLEEVES


POSITION OF THE ANATOMICAL MARKERS

- G corresponds to the axilla (armpit).....
- F corresponds to the widest point of the biceps ...
- E is located at the level of the olecranon
- D equidistant between A and C or at the starting site...... of the œdema
- C to be measured at the narrowest point of the wrist external face in the middle of the radial ulnar joint

Positioning of the measuring tape to measure heights


ARMSLEEVE





The lengths must be measured from the middle of the wrist to the tip of the elbow, then at the outside of the arm up to the desired finished height.

MITTEN






The circumference at the end of the mitten (cA) is located at the head of the metacarpal bones.



Measure the circumference (cX) at the base of the thumb and at the end (cZ). For technical reasons, the length of the thumb (ℓ XZ) cannot be less than 2 cm.



Measure the length of the mitten: 1 / from the wrist (C) to the base of the thumb (X) 2 / from the wrist (C) to the finished height (A)



If possible, please enclose photos of the limb to be fitted. Please draw in the contours of the garment on the diagram and cross unnecessary measures.

FACING VIEW

□ RIGHT ARM □ LEFT ARM

Fill out one form for each side

Models Armsleeve with mitten with thumb (big blocks on the arm and small blocks on the hand) Armsleeve with mitten without thumb (big blocks on the arm and small blocks on the hand)

